

PYRETHROIDS PROFILE

Pesticides in this profile are extracted using a modified QuEChERS (Quick, Easy, Cheap, Effective, Rugged, and Safe) method. Certain analytes may be excluded and some LOQ's may need to be raised due to matrix interferences.

<u>COMPOUND</u>	<u>LOQ (ppm)</u>	<u>COMPOUND</u>	<u>LOQ (ppm)</u>
Acrinathrin	0.01	Permethrin	0.01
Bifenthrin	0.01	Phenothrin	0.01
Cyfluthrin	0.01	Resmethrin	0.01
Cyhalothrin Lambda	0.01	Silafluofen	0.01
Cypermethrin	0.01	Tefluthrin	0.01
Deltamethrin	0.01	Tetramethrin	0.01
Esfenvalerate	0.02	Tralomethrin (screened as Deltamethrin)	0.01
Etofenprox	0.01		
Fenpropathrin	0.01		
Flucythrinate	0.01		
Fluvalinate	0.01		
Halfenprox	0.01		